

# Article The Stream Device—A Retrospective Review of 51 Cases

Karlis Kupcs <sup>1,2</sup>, Patricija Sproge <sup>1</sup>, Katrina Kupca <sup>1</sup> and Pervinder Bhogal <sup>3,\*</sup>

- <sup>1</sup> Department of Neuroradiology, Pauls Stradins Clinical University Hospital, LV-1002 Riga, Latvia; kkupcs@gmail.com (K.K.); patricija.sproge@gmail.com (P.S.); katrina.kupca@gmail.com (K.K.)
- <sup>2</sup> Department of Radiology, Riga Stradins University, LV-1007 Riga, Latvia
- <sup>3</sup> Department of Interventional Neuroradiology, The Royal London Hospital, Barts NHS Trust, London E1 1BB, UK
- \* Correspondence: bhogalweb@aol.com; Tel.: +44-7815937220

Abstract: Mechanical thrombectomy is the gold-standard treatment for patients that have suffered large-vessel occlusion (LVO) stroke. Various different stent-retrievers, aspiration catheters, and techniques have been developed to perform this procedure. We present our initial results regarding the Stream device. Materials and Methods: We performed a retrospective review of a prospectively maintained database at our high-volume centre to identify all patients treated with the Stream device between February 2021 and January 2023. We recorded baseline demographics, NIHSS, ASPECT scores, eTICI scores, complications, and 90-day mRS. Results: We identified 51 patients, 49.0% of whom were male (n = 25), with a median age of 73 (range: 51–89) and a median NIHSS score of 17 (range 4–22), and 68.6% received IV tPA. The median ASPECT score was 10 (range 6–10). Hyperdense clots were seen in 34 cases (66.7%), with a mean clot length of  $12 \pm 6.2$  mm (range 2–26 mm). Clots were located in the anterior circulation in 49 patients. The standard Stream device was used in 78.4% of cases, with Stream 17 being used in 19.6% of cases. The FPE was observed in 25.5% of cases (*n* = 13), with the mFPE being seen in 31.4% of cases (*n* = 16). A final eTICI score of  $\geq$ 2b was achieved in 90.2% of cases (n = 46), and eTICI 2c/3 was seen in 84.3% of cases (n = 43). Furthermore, 24 h CT scans showed that the median ASPECT score was 8 (range 0–10). Good functional outcomes at 90 days (mRS  $\leq$  2) were achieved in 21.6% of cases (n = 11). Conclusions: The Stream device shows acceptable rates of FPE and mFPE compared to existing devices. Further larger studies are required alongside an understanding of the optimal technique for this device's use.

Keywords: stream device; mechanical thrombectomy; ischaemic stroke

## 1. Introduction

Mechanical thrombectomy (MT) has become the gold-standard treatment for acute ischaemic stroke (AIS) caused by large-vessel occlusion (LVO), with studies further demonstrating the efficacy of this procedure in treating patients presenting beyond 6 h [1,2], in relation to posterior circulation, and, more recently, in treating patients with large-core strokes (with ASPECT scores of 3–5) [3–5]. A variety of studies assessing the impact of MT in more distal locations are currently underway, with the expectation that the indications for MT will continue to expand.

Recent advances in our understanding related to clot pathophysiology and structure, clot/stent interactions, and clot/vessel wall interactions have improved our knowledge base regarding ischaemic stroke [6–10]. Over the last decade, there has been a significant improvement in our understanding of the various clot types [6], the interaction between stent-retrievers and clots [7–9], and how both stent-retrievers and clots interact with the vessel wall and in combination during retrieval [10]. This improvement in our understanding has led to the development of stent-retrievers with novel designs [11–15], with the synchronous development of novel aspiration catheters [16–18] and various techniques for optimising the MT procedure and first-pass recanalisation.



Citation: Kupcs, K.; Sproge, P.; Kupca, K.; Bhogal, P. The Stream Device—A Retrospective Review of 51 Cases. J. Clin. Med. 2023, 12, 6384. https://doi.org/10.3390/ jcm12196384

Academic Editors: Benjamin YQ Tan and Franziska Dorn

Received: 11 July 2023 Revised: 17 September 2023 Accepted: 3 October 2023 Published: 6 October 2023



**Copyright:** © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). Manually controlled, expandable, braided devices have been clinically used for the treatment of aneurysms [19–23] and, more recently, cerebral vasospasm [24,25]; however, their development and use in AIS are more recent occurrences. The Stream (Perflow Medical, Tel Aviv, Israel) device is a novel device with a braided structure that, via an actuator at the handle, can be expanded or contracted by the operator.


Currently, there are no publications that have examined the safety and efficacy of this novel device. In this article, we present our initial experience of the use of the Stream device to treat patients with AIS.

## 2. Materials and Methods

# 2.1. The Stream Device

The Stream device is a dynamic neuro-thrombectomy braided net, or an adjustable stent retriever, which is controlled by a proximal handle (Figure 1). At the distal end of the device, there is a 10–13 mm (depending on the device model) atraumatic, radiopaque wire tip. There are 6–8 (depending on the device model) paired wires that form the expandable Cerebral Net<sup>™</sup>, and these can be expanded from 0.4 to 6 mm depending on the model of device chosen. The actuator control handle has two operation modes: auto-lock mode for stepwise radial expansion and free mode for continuous adjustments that can also provide tactile feedback. In auto-lock mode, each 'click' of the handle results in an incremental controlled expansion of the device. There are currently three models:

- Stream;
- Stream XL;
- Stream 17.



**Figure 1.** The Stream device consists of an actuator handle that allows for the precise control of the Cerebral Net ((**A**) and magnified image in (**B**)). At the distal end of the device, there is an atraumatic, radiopaque wire tip (**B**).

Stream and Stream XL are compatible with 0.021-inch-inner-diameter microcatheters, whereas Stream 17 is designed for 0.017-inch-inner-diameter microcatheters. In terms of length, Stream has a maximum braid length of 33 mm, Stream XL has a maximum braid length of 38 mm, and Stream 17 has a maximum braid length of 25 mm. The overall effective length of the Stream is 193 cm.

#### 2.2. Patient Selection

We performed a retrospective review of prospectively maintained databases at a highvolume centre to identify all patients treated with the Stream device between February 2021 and January 2023. The retrospective nature of this study precluded the requirement of ethical approval at the participating centre.

#### 2.3. Study Population

The inclusion criteria included the following:

- Age  $\geq 18$ ;
- National Institutes of Health Stroke Scale (NIHSS) score of ≥4;
- ASPECT score  $\geq$  5;
- LVO on CT angiography;
- Pre-morbid mRS of 0–2;
- Life expectancy of >6 month.

The Stream device was the first device used for treatment of the intracranial occlusion. Patients were excluded if another device or technique (aspiration) was used as the initial strategy for the MT procedure or if there was an incomplete dataset (aside from the 90-day mRS).

All patients underwent non-enhanced CT and CT Angiograms from the arch to the vertex prior to the MT procedure. Patients were eligible for MT if their ASPECT score was  $\geq$ 5, their baseline mRS was  $\leq$ 2, and their life expectancy was greater than 6 months (in the case of a known underlying malignancy), with no upper limit on age. All patients were given IV tPA if they met the criteria.

### 2.4. Endovascular Procedure

Patients were treated either with local anaesthesia or general anaesthesia as per local standard practice.

The exact equipment used varied between the different operators; however, a distal aspiration catheter, typically a 6 Fr Sofia (Microvention, Aliso Viejo, CA, USA), was used in all cases and was typically introduced via a standard guide catheter such as a NeuronMax (Penumbra, Alameda, CA, USA).

Successful reperfusion was defined as eTICI  $\geq$  2b (67%). In addition, first-pass effect (FPE) and modified FPE were defined as eTICI  $\geq$  2c and as eTICI  $\geq$  2b, respectively, after the first thrombectomy attempt.

#### 2.5. Post-Procedure

Post procedure, using non-enhanced CT imaging was routinely performed at  $24 \pm 6$  h unless there was a sudden deterioration in a patient's consciousness level, prompting an earlier scan. In cases of intracranial or extracranial stenting, a CTA was also performed.

The 90-day mRS was recorded via telephone interview or clinic review conducted by a trained stroke physician.

The data recorded included the demographics (age, gender, underlying medical conditions, admission NIHSS score, etc.), use of IV tPA, radiological findings (including the ASPECT score prior to MT), clot location, endovascular procedural information, and relevant timing. The eTICI score was recorded after the first pass of the device and at the end of the procedure alongside any complications and the 90-day mRS where available.

#### 3. Results

#### 3.1. Baseline Demographics

In total, 51 patients met our inclusion and exclusion criteria. The median age of the participants was 73 (range 51–89), and 49.0% were male (n = 25). Pre-existing hypertension was common (n = 40, 78.4%), diabetes mellitus was less common (n = 13, 25.5%), and atrial fibrillation was noted in just over half of the patients (n = 26, 51%). All patients were recorded as having a baseline mRS of 0 prior to the MT procedure. The results are summarised in Table 1.

Baseline Data	n = 51
Demographics	
Age	Median 73 (range 51–89)
Female	51% ( <i>n</i> = 26)
Co-Morbidities	
Smoking	20 (39.2%)
Hypertension	40 (78.4%)
Diabetes Mellitus	13 (25.5%)
Atrial Fibrillation	26 (51%)
Pre-Morbid mRS	
0	51 (100%)

Table 1. Baseline demographics and medical data.

## 3.2. Clinical and Angiographic Results

The median NIHSS score at presentation was 17 (range 4–22), and 68.6% received IV tPA prior to MT. The most common suspected cause of a stroke was thought to be cardioembolic (52.9%). Right-sided LVOs were more frequently seen (56.9%), and eight cases involved tandem lesions (15.7%). The median ASPECT score upon conducting a plain CT scan was 10 (range 6–10). Hyperdense clots were seen in 34 cases (66.7%), with a mean clot length of  $12 \pm 6.2$  mm (range 2–26 mm). The vast majority of clots were located in the anterior circulation, with only two posterior basilar occlusions (3.9%).


The results are summarized in Table 2.

**Table 2.** Baseline clinical and imaging data, including suspected cause of stroke, NIHSS score, clot location, and ASPECT score.


Stroke Data	
NIHSS	Median 17 (range 4–22)
IV tPA	35 (68.6%)
Suspected Cause	
Cardioembolic	27 (52.9%)
Large-artery Atherosclerosis	6 (11.8%)
Mixed	5 (9.8%)
ESUS (Embolic Stroke of Undetermined Source)	13 (25.5%)
Imaging Findings	
Side	
R	29 (56.9%)
L	20 (39.2%)
Mid	2 (3.9%)
Tandem Lesion	
Y	8 (15.7%)
Clot Location	
ICA	16 (31.4%)
M1	29 (56.9%)
M2	4 (7.8%)
BA	2 (3.9%)
Hyperdense Clot	34 (66.7%)
Clot Length	$12 \pm 6.2$ mm (range: 2–26 mm)
ASPECT Score	Median 10 (range: 6–10)

## 3.3. Procedural Outcomes

The vast majority of patients were operated on under local anaesthesia (n = 49, 96.1%). The standard Stream device was used as the first-line device in the majority of procedures (n = 40, 78.4%), with the smaller Stream 17 used in 19.6% of cases (n = 10). The first-pass effect (FPE) was seen in 25.5% of cases (n = 13), with a modified FPE seen in 31.4% (n = 16) (Figures 2 and 3). Distal embolisation was seen in 5.9% (n = 3) of cases, and an embolus in new territory was not seen in any cases.



**Figure 2.** A male patient with 2 h history of acute left-sided weakness and an NIHSS score of 17. Non-contrast CT scan (**A**) revealed some subtle loss of grey–white matter differentiation involving the lentiform nucleus (ASPECT score 9), and CT angiography confirmed a right terminal ICA and M1 occlusion (**B**).



**Figure 3.** Angiography confirmed the terminal right ICA and M1 occlusion (**A**). The thrombus was crossed using a Headway 21 microcatheter, and a Stream device was expanded for 5 min (**B**). Withdrawal into the aspiration catheter was performed, during which the first-pass effect was observed (eTICI 2c) (**C**). Non-contrast CT scan of the head performed at 24 h demonstrated no evidence of haemorrhage and an ASPECT score of 8 (**D**).

In 19.6% (n = 10) of cases, Stream could not be used to remove the occlusion, and bailout with alternative devices/techniques was required. An intracranial stent was permanently implanted in one of these cases (1.9%). Carotid stents were implanted in two cases (3.9%).

The median number of passes when only Stream was used was 2 (range: 1-4), with a mean number of passes equal to  $2.0 \pm 0.92$ . At the end of the procedure, after all the devices and techniques had been used, a final eTICI score of  $\geq$ 2b was achieved in 90.2% of cases (n = 46), with eTICI 2c/3 seen in 84.3% (n = 43) of cases. The results are summarised in Table 3.

**Table 3.** Procedural data including the model of Stream used, type of anaesthesia, angiographic result following MT, and distal and new-territory embolisation.

Procedural Data	
Anaesthesia	
LA	49 (96.1%)
GA	2 (3.9%)
Distal Aspiration Catheter	51 (100%)
Stream	
Stream	40 (78.4%)
Stream XL	1 (2%)
Stream 17	10 (19.6%)
Angiographic Results	
Median No. of passes	2 (range 1–8)
Median No. of passes with Stream	2 (range 1–4)
FPE (eTICI 2c/3)	13 (25.5%)
Modified FPE (eTICI $\geq 2b$ )	16 (31.4%)
Final eTICI	
0–2a	5 (9.8%)
2b	3 (5.9%)
2c	4 (7.8%)
3	39 (76.5%)
Distal Embolisation	3 (5.9%)
Embolisation to New territory	0
Bailout Required	<i>n</i> = 10
Resistant Clot	5 (9.8%)
Device Damage	4 (7.8%)
Intracranial Stent Implanted	1 (1.9%)
Carotid Stent Implanted	2 (3.9%)

## 3.4. Follow-up Imaging and Clinical Results

A follow-up CT scan performed at 24 h was available for all patients. The median ASPECT score was 8 (range 0–10). Subarachnoid haemorrhage was seen in eight patients (15.7%), and seven patients developed sICH (13.7%), three of whom also had SAH. The 90-day mRS was available for all patients, and a good functional outcome (mRS 0–2) was achieved in 21.6% of cases (n = 11), with an mRS of 6 in 43.1% of cases. The results are summarised in Table 4.

Table 4. Angiographic and clinical follow-up data.

Median 8 (range 0–10)
7(1270)
7 (13.7%)
8 (15.7%)
11 (21.6%)
18 (35.3%)
22 (43.1%)

#### 4. Discussion

In our initial experience with the Stream device, we observed a similar rate of efficacy with respect to FPE and mFPE when compared to other frequently used stent-retrievers such as the Solitaire (Medtronic) [26–29] and Trevo (Stryker) stent-retrievers [30–32], whose efficacy averages between 25 and 40%.

The Stream device has several potential technical advantages that come into effect during a mechanical thrombectomy. Unlike standard stent-retrievers, this device is under the direct control of the operating physician in the sense that it does not automatically expand like a standard device. Therefore, if the Stream device has been unsheathed in a suboptimal position, it can be repositioned to optimise the clot/device interaction. Although one can achieve the same result with standard stent-retrievers, a potential drawback is that they will immediately expand upon unsheathing, and this will not only compress the clot but also result in microthrombi that have the potential to travel distally [10]. Similarly, the closed nature of the distal end of the Stream device's braided net may also act as a distal clot catcher that prevents, at least to a degree, the loss of distal thrombi during the expansion of the device, as was recently demonstrated with another device [33]. Interestingly, in our study, distal emboli were only seen in 5.9% of cases despite having been reported in  $\approx 25\%$  of cases for devices with distal clot catchers [12]. It is also noteworthy that the average clot length in the three cases where distal clot embolisation occurred was  $17.7\pm6.8$  mm, which was considerably longer than the average clot length of  $11.6\pm6.1$  mm in the cases in which no distal embolisation was seen. This may be due to the fact that distal embolisation may, at least in part, be due to the stent-retriever's length-to-clot-length ratio as has been suggested by Belachew et al. [34]. The fact that the Stream device shortens during expansion may mean that for longer clots, a longer device should be chosen in order to optimise the device/clot interaction, and further studies are required in order to determine if this holds true for the Stream device as it appears to for other devices [35–37].

The radial force of a stent-retriever is generally fixed and decreases at increasing diameters. Although techniques for increasing the radial force of standard stent-retrievers have been developed, such as the 'push and fluff' technique, these are unlikely to lead to very significant increases in the radial force developed. It has previously been shown that the migration of stent-struts into a clot is determined, at least in part, by the radial force exerted as well as the size of the pores [38]. The radially expandable devices do not share this limitation, and, in fact, the radial force can be increased with manual expansion in a similar manner to, but to a much lesser degree, angioplasty balloons. Conversely, the radial force can also be decreased during the procedure in case excessive resistance is felt. This allows the operator to expand and size the device according to the target vessel; importantly, this can be adapted during the thrombus's retrieval. This feature may be of particular importance in platelet-rich white thrombi, which prove resistant to standard mechanical thrombectomy approaches. Even in the event of a failure to remove these clots, valuable information can be gleaned from the attempt, as the Stream device is visible along its entire length and, as such, the exact location of the clot can be determined if bailout angioplasty and stenting are to be considered.

In our experience, we have only performed a standard unsheathing and expansion of the Stream device; however, there is some evidence to suggest that repeated inflation and deflation may offer an advantage. In the study conducted by Kara et al. [39], using the Tigertriever, the first-pass recanalisation (mTICI  $\geq 2b$ ) rate using a repeated inflation and deflation technique was 47.8% compared to 31.6% when using a standard unsheathing and expansion technique. The authors believe that this might have been due to a combination of factors. Repeated inflation and over inflation may result in the greater enlargement of the cells, allowing more clot fragments to be entrapped within the lumen of the device. The authors also suggest that oversizing causes better apposition of the stent to the vessel wall, which enhances stent–clot interaction. We believe that a further potential reason for this improved recanalisation can be the fact that static expansion will not only potentially compress the clot but also increase the frictional forces between the clot and the vessel

wall, whereas repeated inflation and deflation may 'pull' the clot away from the vessel wall with each deflation and hence aid in overcoming these frictional forces during retrieval. It remains to be seen whether the repeated expansion technique results in improved recanalisation rates; however, this warrants further investigation with larger cohorts as well as bench-side studies.

Our study has several limitations that are inherent to retrospective observational studies. We did not follow a predefined study protocol that included a fixed number of attempts with the Stream device before switching to a different device/technique, and we did not standardise the other pieces of equipment used. The eTICI scores and follow-up imaging procedures were not adjudicated by an independent core lab.

#### 5. Conclusions

The Stream device showed acceptable rates of FPE and mFPE compared to existing devices. Further larger studies are required alongside an understanding of the optimal technique for this device's use.

Author Contributions: K.K. (Karlis Kupcs)—data collection, editing, and guarantor; P.S.—data collection, editing, and revision; K.K. (Katrina Kupca)—data collection, editing, and revision; P.B.—manuscript preparation, editing, and revision. All authors have read and agreed to the published version of the manuscript.

Funding: This research received no external funding.

**Institutional Review Board Statement:** Ethics committee approval is not required for a retrospective study of this nature.

**Informed Consent Statement:** Informed consent was obtained from all individual patients who were able to provide consent. For patients that could not provide consent, the thombectomy procedures were discussed with their next of kin/available family members and consent was obtained by the operating team and stroke physicians.

**Conflicts of Interest:** P.B.—consulting agreements Phenox, Cerenovus, Balt, Neurovasc Technologies, Brainomix, Perflow Medical, Perfuze, PockIt Diagnositcs, Vesalio, and BT. The remaining authors declare no conflict of interest.

## References

- Nogueira, R.G.; Jadhav, A.P.; Haussen, D.C.; Bonafe, A.; Budzik, R.F.; Bhuva, P.; Yavagal, D.R.; Ribo, M.; Cognard, C.; Hanel, R.A.; et al. Thrombectomy 6 to 24 Hours after Stroke with a Mismatch between Deficit and Infarct. *N. Engl. J. Med.* 2018, 378, 11–21. [CrossRef] [PubMed]
- Albers, G.W.; Marks, M.P.; Kemp, S.; Christensen, S.; Tsai, J.P.; Ortega-Gutierrez, S.; McTaggart, R.A.; Torbey, M.T.; Kim-Tenser, M.; Leslie-Mazwi, T.; et al. Thrombectomy for Stroke at 6 to 16 Hours with Selection by Perfusion Imaging. N. Engl. J. Med. 2018, 378, 708–718. [CrossRef] [PubMed]
- Sarraj, A.; Hassan, A.E.; Abraham, M.G.; Ortega-Gutierrez, S.; Kasner, S.E.; Hussain, M.S.; Chen, M.; Blackburn, S.; Sitton, C.W.; Churilov, L.; et al. Trial of Endovascular Thrombectomy for Large Ischemic Strokes. *N. Engl. J. Med.* 2023, 388, 1259–1271. [CrossRef] [PubMed]
- Uchida, K.; Shindo, S.; Yoshimura, S.; Toyoda, K.; Sakai, N.; Yamagami, H.; Matsumaru, Y.; Matsumoto, Y.; Kimura, K.; Ishikura, R.; et al. Association Between Alberta Stroke Program Early Computed Tomography Score and Efficacy and Safety Outcomes With Endovascular Therapy in Patients With Stroke From Large-Vessel Occlusion: A Secondary Analysis of the Recovery by Endovascular Salvage for Cerebral Ultra-Acute Embolism-Japan Large Ischemic Core Trial (RESCUE-Japan LIMIT). *JAMA Neurol.* 2022, 79, 1260–1266. [CrossRef] [PubMed]
- Huo, X.; Ma, G.; Tong, X.; Zhang, X.; Pan, Y.; Nguyen, T.N.; Yuan, G.; Han, H.; Chen, W.; Wei, M.; et al. Trial of Endovascular Therapy for Acute Ischemic Stroke with Large Infarct. *N. Engl. J. Med.* 2023, 388, 1272–1283. [CrossRef] [PubMed]
- Staessens, S.; Denorme, F.; Francois, O.; Desender, L.; Dewaele, T.; Vanacker, P.; Deckmyn, H.; Vanhoorelbeke, K.; Andersson, T.; De Meyer, S.F. Structural Analysis of Ischemic Stroke Thrombi: Histological Indications for Therapy Resistance. *Haematologica* 2020, 105, 498–507. [CrossRef] [PubMed]
- Dumitriu LaGrange, D.; Bernava, G.; Reymond, P.; Wanke, I.; Vargas, M.I.; Machi, P.; Lövblad, K.-O. A High Resolution Scanning Electron Microscopy Analysis of Intracranial Thrombi Embedded along the Stent Retrievers. *Sci. Rep.* 2022, 12, 8027. [CrossRef] [PubMed]

- Tsumoto, T.; Tsurusaki, Y.; Tokunaga, S. Interaction between the Stent Strut and Thrombus Characterized by Contrast-Enhanced High-Resolution Cone Beam CT during Deployment of the Solitaire Stent Retriever. J. Neurointerv. Surg. 2017, 9, 843–848.
   [CrossRef]
- Autar, A.S.A.; Hund, H.M.; Ramlal, S.A.; Hansen, D.; Lycklama, À.; Nijeholt, G.J.; Emmer, B.J.; de Maat, M.P.M.; Dippel, D.W.J.; van der Lugt, A.; et al. High-Resolution Imaging of Interaction Between Thrombus and Stent-Retriever in Patients with Acute Ischemic Stroke. J. Am. Heart Assoc. 2018, 7, e008563. [CrossRef]
- Liu, Y.; Gebrezgiabhier, D.; Reddy, A.S.; Davis, E.; Zheng, Y.; Arturo Larco, J.L.; Shih, A.J.; Pandey, A.S.; Savastano, L.E. Failure Modes and Effects Analysis of Mechanical Thrombectomy for Stroke Discovered in Human Brains. *J. Neurosurg.* 2021, 136, 197–204. [CrossRef]
- Brouwer, P.A.; Yeo, L.L.L.; Holmberg, A.; Andersson, T.; Kolloch, J.; KuntzeSöderqvist, Å.; Ohlsson, M.; Holmin, S.; Mpotsaris, A.; Gontu, V.K.; et al. Thrombectomy Using the EmboTrap Device: Core Laboratory-Assessed Results in 201 Consecutive Patients in a Real-World Setting. *J. Neurointerv. Surg.* 2018, 10, 964–968. [CrossRef] [PubMed]
- Omrani, O.; Hafezi-Bakhtiari, N.; DeSouza, P.; Nikola, C.; Wong, K.; Lansley, J.; Dhillon, P.; Makalanda, L.; Chan, N.; Harrison, T.; et al. The Initial Experience with the Embotrap III Stent-Retriever in a Real World Setting. *Interv. Neuroradiol. J. Peritherapeutic Neuroradiol. Surg. Proced. Relat. Neurosci.* 2022, 15910199221142096. [CrossRef] [PubMed]
- Zaidat, O.O.; Bozorgchami, H.; Ribó, M.; Saver, J.L.; Mattle, H.P.; Chapot, R.; Narata, A.P.; Francois, O.; Jadhav, A.P.; Grossberg, J.A.; et al. Primary Results of the Multicenter ARISE II Study (Analysis of Revascularization in Ischemic Stroke With EmboTrap). *Stroke* 2018, 49, 1107–1115. [CrossRef] [PubMed]
- Bajrami, A.; Ertugrul, O.; Senadim, S.; Erdem, E.; Baltacioglu, F.; Geyik, S. First Pass Results of Mechanical Thrombectomy with Two-Drop Zone NeVaTM Device. *Interv. Neuroradiol. J. Peritherapeutic Neuroradiol. Surg. Proced. Relat. Neurosci.* 2022, 15910199221135308. [CrossRef] [PubMed]
- Akpinar, C.K.; Ozdemir, A.O.; Gurkas, E.; Bilgic, A.B.; Aykac, O.; Inanc, Y.; Giray, S. Favorable First-Pass Recanalization Rates with NeVa<sup>TM</sup> Thrombectomy Device in Acute Stroke Patients: Initial Clinical Experience. *Interv. Neuroradiol. J. Peritherapeutic Neuroradiol. Surg. Proced. Relat. Neurosci.* 2021, 27, 107–113. [CrossRef]
- 16. Makalanda, L.; Lansley, J.; Wong, K.; Spooner, O.; Bhogal, P. The Q and A-The MIVI Q Catheters for Aspiration Thrombectomy-Initial Experience from London. *J. Clin. Med.* **2021**, *10*, 5844. [CrossRef] [PubMed]
- de Rueda, M.E.; Ballenilla Marco, F.; Garmendia Lopetegui, E.; Pumar, J.M.; Zamarro, J.; García-Villalba, B.; Díaz-Pérez, J.; Mosqueira, A.; Lüttich, A.; Larrea, J.-A.; et al. Thrombectomy Aspiration Post-Market Study in Acute Stroke with the Q Aspiration Catheter: The TAPAS Study. J. Neurointerv. Surg. 2022, 15, 674–678. [CrossRef] [PubMed]
- Kobeissi, H.; Ghozy, S.; Flood, R.; Mortimer, A.; Crossley, R.; Cox, A.; Minks, D.; Wareham, J. Mechanical Thrombectomy with Q Catheter in Stroke Caused by Primary and Secondary Distal and Medium Vessel Occlusions. *Interv. Neuroradiol. J. Peritherapeutic Neuroradiol. Surg. Proced. Relat. Neurosci.* 2023, 15910199231167916. [CrossRef]
- Sirakov, S.; Sirakov, A.; Minkin, K.; Karakostov, V.; Raychev, R. Early Clinical Experience with Cascade: A Novel Temporary Neck Bridging Device for Embolization of Intracranial Aneurysms. J. Neurointerv. Surg. 2020, 12, 303–307. [CrossRef]
- Tomasello, A.; Hernandez, D.; Gramegna, L.L.; Aixut, S.; Barranco Pons, R.; Jansen, O.; Zawadzki, M.; Lopez-Rueda, A.; Parra-Fariñas, C.; Piñana, C.; et al. Early Experience with a Novel Net Temporary Bridging Device (Cascade) to Assist Endovascular Coil Embolization of Intracranial Aneurysms. J. Neurosurg. 2020, 134, 591–599. [CrossRef]
- 21. Bhogal, P.; Wong, K.; Makalanda, H. The Cascade Device—In Vitro Tests to Assess Coil Protrusion into the Parent Vessel. Interv. Neuroradiol. J. Peritherapeutic Neuroradiol. Surg. Proced. Relat. Neurosci. 2020, 26, 494–500. [CrossRef] [PubMed]
- 22. Fischer, S.; Weber, A.; Carolus, A.; Drescher, F.; Götz, F.; Weber, W. Coiling of Wide-Necked Carotid Artery Aneurysms Assisted by
- a Temporary Bridging Device (Comaneci): Preliminary Experience. *J. Neurointerv. Surg.* 2017, *9*, 1039–1097. [CrossRef] [PubMed]
  23. Gupta, R.; Kolodgie, F.D.; Virmani, R.; Eckhouse, R. Comaneci Neck Bridging Device for the Treatment of Cerebral Aneurysms. *J. Neurointerv. Surg.* 2016, *8*, 181–185. [CrossRef] [PubMed]
- Bhogal, P.; Simpanen, T.; Wong, K.; Bushi, D.; Sirakov, M.A.; Sirakov, S.; Aggour, M.; Makalanda, L. Use of the Cascade Expandable Net to Treat Cerebral Vasospasm—Initial Clinical Experience from a Single Centre with in Vitro Benchside Tests. *CVIR Endovasc.* 2021, 4, 82. [CrossRef] [PubMed]
- Salem, M.M.; Khalife, J.; Desai, S.; Sharashidze, V.; Badger, C.; Kuhn, A.L.; Monteiro, A.; Salahuddin, H.; Siddiqui, A.H.; Singh, J.; et al. COManeci MechANical Dilation for Vasospasm (COMMAND): Multicenter Experience. *J. Neurointerv. Surg.* 2022, 15, 864–870. [CrossRef] [PubMed]
- Zaidat, O.O.; Castonguay, A.C.; Gupta, R.; Sun, C.-H.J.; Martin, C.; Holloway, W.E.; Mueller-Kronast, N.; English, J.D.; Linfante, I.; Dabus, G.; et al. North American Solitaire Stent Retriever Acute Stroke Registry: Post-Marketing Revascularization and Clinical Outcome Results. J. Neurointerv. Surg. 2014, 6, 584–588. [CrossRef] [PubMed]
- Dávalos, A.; Pereira, V.M.; Chapot, R.; Bonafé, A.; Andersson, T.; Gralla, J. Solitaire Group Retrospective Multicenter Study of Solitaire FR for Revascularization in the Treatment of Acute Ischemic Stroke. *Stroke J. Cereb. Circ.* 2012, 43, 2699–2705. [CrossRef] [PubMed]
- 28. Broussalis, E.; Trinka, E.; Hitzl, W.; Wallner, A.; Chroust, V.; Killer-Oberpfalzer, M. Comparison of Stent-Retriever Devices versus the Merci Retriever for Endovascular Treatment of Acute Stroke. *AJNR Am. J. Neuroradiol.* **2013**, *34*, 366–372. [CrossRef]

- Campbell, B.C.V.; Hill, M.D.; Rubiera, M.; Menon, B.K.; Demchuk, A.; Donnan, G.A.; Roy, D.; Thornton, J.; Dorado, L.; Bonafe, A.; et al. Safety and Efficacy of Solitaire Stent Thrombectomy: Individual Patient Data Meta-Analysis of Randomized Trials. *Stroke* 2016, 47, 798–806. [CrossRef]
- Binning, M.J.; Bartolini, B.; Baxter, B.; Budzik, R.; English, J.; Gupta, R.; Hedayat, H.; Krajina, A.; Liebeskind, D.; Nogueira, R.G.; et al. Trevo 2000: Results of a Large Real-World Registry for Stent Retriever for Acute Ischemic Stroke. *J. Am. Heart Assoc.* 2018, 7, e010867. [CrossRef]
- Jadhav, A.P.; Desai, S.M.; Budzik, R.F.; Gupta, R.; Baxter, B.; English, J.D.; Bartolini, B.M.; Krajina, A.; Haussen, D.C.; Nogueira, R.G.; et al. First Pass Effect in Patients with Large Vessel Occlusion Strokes Undergoing Neurothrombectomy: Insights from the Trevo Retriever Registry. J. Neurointerv. Surg. 2021, 13, 619–622. [CrossRef] [PubMed]
- 32. Protto, S.; Pienimäki, J.-P.; Seppänen, J.; Matkaselkä, I.; Ollikainen, J.; Numminen, H.; Sillanpää, N. TREVO and Capture LP Have Equal Technical Success Rates in Mechanical Thrombectomy of Proximal and Distal Anterior Circulation Occlusions. *J. Neurointerv. Surg.* **2017**, *9*, 644–649. [CrossRef] [PubMed]
- 33. Li, J.; Tiberi, R.; Bhogal, P.; Buhk, J.-H.; Behme, D.; Tomasello, A.; Ribo, M. Impact of Stent-Retriever Tip Design on Distal Embolization during Mechanical Thrombectomy: A Randomized in Vitro Evaluation. *J. Neurointerv. Surg.* **2023**. [CrossRef]
- Belachew, N.F.; Dobrocky, T.; Meinel, T.R.; Hakim, A.; Vynckier, J.; Arnold, M.; Seiffge, D.J.; Wiest, R.; Piechowiak, E.I.; Fischer, U.; et al. Risks of Undersizing Stent Retriever Length Relative to Thrombus Length in Patients with Acute Ischemic Stroke. *AJNR Am.* J. Neuroradiol. 2021, 42, 2181–2187. [CrossRef] [PubMed]
- 35. Serna Candel, C.; Aguilar Pérez, M.; Bäzner, H.; Henkes, H.; Hellstern, V. First-Pass Reperfusion by Mechanical Thrombectomy in Acute M1 Occlusion: The Size of Retriever Matters. *Front. Neurol.* **2021**, *12*, 679402. [CrossRef] [PubMed]
- Haussen, D.C.; Al-Bayati, A.R.; Grossberg, J.A.; Bouslama, M.; Barreira, C.; Bianchi, N.; Frankel, M.R.; Nogueira, R.G. Longer Stent Retrievers Enhance Thrombectomy Performance in Acute Stroke. J. Neurointerv. Surg. 2019, 11, 6–8. [CrossRef] [PubMed]
- Zaidat, O.O.; Haussen, D.C.; Hassan, A.E.; Jadhav, A.P.; Mehta, B.P.; Mokin, M.; Mueller-Kronast, N.H.; Froehler, M.T. Impact of Stent Retriever Size on Clinical and Angiographic Outcomes in the STRATIS Stroke Thrombectomy Registry. *Stroke* 2019, 50, 441–447. [CrossRef]
- Wenger, K.; Nagl, F.; Wagner, M.; Berkefeld, J. Improvement of Stent Retriever Design and Efficacy of Mechanical Thrombectomy in a Flow Model. *Cardiovasc. Intervent. Radiol.* 2013, 36, 192–197. [CrossRef]
- Kara, B.; Selcuk, H.H.; Erbahceci Salik, A.; Zalov, H.; Yildiz, O.; Gul, G.; Balkan, B. Single-Center Experience with the Tigertriever Device for the Recanalization of Large Vessel Occlusions in Acute Ischemic Stroke. J. Neurointerv. Surg. 2019, 11, 455–459. [CrossRef]

**Disclaimer/Publisher's Note:** The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.